EQUIPMENT:

MULTI-PURPOSE WORKSTATION™ STABILIZER

MULTI-PURPOSE WORKSTATION™ STABILIZER LP

Symbol	Title of Symbol	Description of Symbol
	Manufacturer (ISO 15223-1, 5.1.1)	Indicates the medical device manufacturer.
	Authorized Representative in the European Community (ISO 15223-1, 5.1.2)	Indicates the Authorized Representative in the European Community.
	Date of manufacture (ISO 15223-1, 5.1.3)	Indicates the date when the medical device was manufactured.
	Batch code (ISO 15223-1, 5.1.5)	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog number (ISO 15223-1, 5.1.6)	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	Serial number (ISO 15223-1, 5.1.7)	Indicates the manufacturer's serial number so that a specific medical device can be identified.
	Consult instructions for use (ISO 15223-1, 5.4.3)	Indicates the need for the user to consult the instructions for use.
	European Conformity (EU MDR 2017/745, Article 20)	Indicates manufacturer declaration that the product complies with the essential requirements of the relevant European health, safety and environment protection legislation.
	Medical Device (MedTech Europe Guidance: Use of Symbols to Indicate Compliance with the MDR)	Indicates the product is a medical device.
	Quantity (IEC 60878, 2794)	To indicate the number of pieces in the package.

GENERAL EQUIPMENT INFORMATION

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

WARNING

- Before use, you should be trained in ultrasonography. For instructions on the use of your transducer, see your system's user guide.
- Prior to use, read and understand all instructions and warnings.
- Prior to use inspect device for signs of damage, if damage is evident do not use.
- Equipment is not to be modified without CIVCO authorization.
- Stabilizer is designed and validated for use with CIVCO accessories. For CIVCO Steppers and other accessories, visit www.CIVCO.com.
- Stabilizer is packaged non-sterile and is reusable. To avoid possible patient contamination, ensure stabilizer is properly cleaned and disinfected before each use. See reprocessing section for instructions on how to properly clean and disinfect.
- For illustration purposes only, equipment may be shown without a drape. Always place a drape over the equipment to protect patients and users from cross-contamination.
- If the product malfunctions during use or is no longer able to achieve its intended use, stop using the product and contact CIVCO.
- Report serious incidents related to the product to CIVCO and the competent authority in your Member State, or appropriate regulatory authorities.

INTENDED USE

The equipment is intended to provide fixation, support and manipulation of transrectal ultrasound imaging probes during insertion and final placement.

INDICATIONS FOR USE

Prostate - Diagnostic imaging and minimally invasive puncture procedures.

Surgical (Prostate) - Diagnostic imaging and puncture procedures.

PATIENT POPULATION

Equipment is for use in adult males with suspected or diagnosed prostate cancer.

INTENDED USERS

Equipment should be used by clinicians medically trained in ultrasound imaging. User groups may include, but are not limited to: Physicists, Radiation Oncologists, Surgeons and Urologists.


PERFORMANCE CHARACTERISTICS

Stabilizer features single-point locking mechanism to rapidly lock the stabilizer position without transducer migration.


NOTE: For a summary of clinical benefits for this product, visit www.CIVCO.com.

MULTI-PURPOSE WORKSTATION™ STABILIZER

ASSEMBLING FLOORSTAND

USING FLOORSTAND

- I. Vertical Adjustment Knob
- II. Arm Locking Knob
- III. Hex Shaft Attachment
- IV. Vertical Adjustment Shaft
- V. Wheel Locking Lever

1. Unlock all wheel locking levers to move workstation floorstand into optimal working position.
2. When optimal working position, push wheel locking levers down to secure.
3. To adjust optimal working height for horizontal swivel joint, first loosen vertical adjustment knob then lift vertical adjustment shaft to desired height. Tighten vertical adjustment knob to secure.

 **WARNING**


- Ensure mounting system is secure.

4. Attach desired accessory to hex shaft.
5. Loosen arm locking knob to position arm.
6. Tighten arm locking knob to secure position.

MULTI-PURPOSE WORKSTATION™ STABILIZER LP (SINGLE-SIDED MOUNT)



SETTING UP ARM (SINGLE-SIDED MOUNT)

- Loosen rail clamp levers on table mount. Slide table mount onto table rail and tighten rail clamps levers to secure. Center Y adjustment.

- I. Stabilizer Quick Connect Knob
- II. Knob Y
- III. Rail Clamp
- IV. Table Rail
- V. Rail Clamp Levers
- VI. Stabilizer Mounting Plate
- VII. Hex Shaft Attachment
- VIII. Arm Locking Knob

- Fit dovetailed mounting plate to mating stabilizer mounting plate on table mount. Secure with stabilizer quick connect knob.

FUNCTIONAL PREUSE CHECKS

PRIOR TO USE

- Attach table mount at a suitable clamping point on table rails to ensure it remains securely in position.
- Ensure stabilizer moves freely when locking knob is loosened and secures firmly when locking knob is tightened.
- Perform table clearance check. Before indexing the table, make sure stabilizer and table mount do not interfere.

REPROCESSING

WARNING

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.
- See your system's user guide for reprocessing transducer between uses.
- These procedures for reprocessing have been validated for effectiveness and compatibility. Equipment could be damaged or cross-contaminated due to improper reprocessing.
- Do not place in a mechanical washer.
- Do not use ethylene oxide or autoclave to sterilize.

- Disinfect all surfaces with common germicidal or antiseptic solution, such as alcohol or hydrogen peroxide.
- Wipe dry with lint-free cloth or allow to air dry prior to reassembly and use.

STORAGE AND DISPOSAL

WARNING

- Dispose of single-use components as infectious waste. Reprocess reusable components after each use according to instructions for use.

NOTE:

- Device should be stored to keep all components together and secure.
- For questions or to order additional CIVCO products, please call +1 319-248-6757 or 1-800-445-6741 or visit www.CIVCO.com.
- All product to be returned must be in its original packaging. Contact CIVCO for further instructions as needed.

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

STABILIZATOR MULTI-PURPOSE WORKSTATION™

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


**Australian
Sponsor Address**

SCANMEDICS
Unit 6, 15-17 Gibbes St.
St. Chatswood, NSW 2067
Australia

CIVCO Medical Instruments Co., Inc.
102 First Street South, Kalona, IA 52247-9589
319.248.6757, 800.445.6741
Fax: 319.248.6660, 877.329.2482

COPYRIGHT © 2020 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ACCUCARE AND MULTI-PURPOSE WORKSTATION ARE TRADEMARKS OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS.
PRINTED IN USA. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW.

WWW.CIVCO.COM

2020-08-17 043-659-5A

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MPS GmbH