
EQUIPMENT:

MICRO-TOUCH® STABILIZER (DUAL-SIDED TABLE MOUNT)

MICRO-TOUCH® STABILIZER LP (SINGLE-SIDED TABLE MOUNT)

Symbol	Title of Symbol	Description of Symbol
	Manufacturer (ISO 15223-1, 5.1.1)	Indicates the medical device manufacturer.
	Authorized Representative in the European Community (ISO 15223-1, 5.1.2)	Indicates the Authorized Representative in the European Community.
	Date of manufacture (ISO 15223-1, 5.1.3)	Indicates the date when the medical device was manufactured.
	Batch code (ISO 15223-1, 5.1.5)	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog number (ISO 15223-1, 5.1.6)	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	Serial number (ISO 15223-1, 5.1.7)	Indicates the manufacturer's serial number so that a specific medical device can be identified.
	Consult instructions for use (ISO 15223-1, 5.4.3)	Indicates the need for the user to consult the instructions for use.
	European Conformity (EU MDR 2017/745, Article 20)	Indicates manufacturer declaration that the product complies with the essential requirements of the relevant European health, safety and environment protection legislation.
	Medical Device (MedTech Europe Guidance: Use of Symbols to Indicate Compliance with the MDR)	Indicates the product is a medical device.
	Quantity (IEC 60878, 2794)	To indicate the number of pieces in the package.

GENERAL EQUIPMENT INFORMATION

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

WARNING

- Before use, you should be trained in ultrasonography. For instructions on the use of your transducer, see your system's user guide.
- Prior to use, read and understand all instructions and warnings.
- Prior to use inspect device for signs of damage, if damage is evident do not use.
- Equipment is not to be modified without CIVCO authorization.
- Stabilizer is designed and validated for use with CIVCO accessories. For CIVCO Steppers and other accessories, visit www.CIVCO.com.
- Stabilizer is packaged non-sterile and is reusable. To avoid possible patient contamination, ensure stabilizer is properly cleaned and disinfected before each use. See reprocessing section for instructions on how to properly clean and disinfect.
- For illustration purposes only, equipment may be shown without a drape. Always place a drape over the equipment to protect patients and users from cross-contamination.
- If the product malfunctions during use or is no longer able to achieve its intended use, stop using the product and contact CIVCO.
- Report serious incidents related to the product to CIVCO and the competent authority in your Member State, or appropriate regulatory authorities.

INTENDED USE

The equipment is intended to provide fixation, support and manipulation of transrectal ultrasound imaging probes during insertion and final placement.

INDICATIONS FOR USE

Prostate - Diagnostic imaging and minimally invasive puncture procedures.

Surgical (Prostate) - Diagnostic imaging and puncture procedures.

PATIENT POPULATION

Equipment is for use in adult males with suspected or diagnosed prostate cancer.

INTENDED USERS

Equipment should be used by clinicians medically trained in ultrasound imaging. User groups may include, but are not limited to: Physicists, Radiation Oncologists, Surgeons and Urologists.

PERFORMANCE CHARACTERISTICS


Stabilizer features stable single-arm fixation with a single-point locking mechanism to instantly fix the stabilizer in the desired position.

NOTE: For a summary of clinical benefits for this product, visit www.CIVCO.com.

ASSEMBLING TRANSPORTATION STAND


For use with Micro-Touch® Stabilizer (Dual-Sided Table Mount)

1. Insert threaded stud end of tube through center of rolling base.
2. Slide locating washer over threaded stud and nest into center of hole in rolling base.
3. Turn locking knob onto threaded stud and secure firmly.

- I. Threaded stud
- II. Rolling base
- III. Locating washer
- IV. Locking knob


4. The transportation stand provides a convenient and mobile platform to allow ease of cleaning and handling within the surgical environment. Stand shown with stabilizer and table mount clamped in place.

- I. Stabilizer
- II. Table mount
- III. Transportation stand

SETTING UP MICRO-TOUCH® STABILIZER (DUAL-SIDED TABLE MOUNT)

1. Loosen all knobs and levers on table mount before removing from transportation stand. Grasp rail clamps and slide entire apparatus into position on end of table. Lock rail clamps onto table by tightening rail clamp levers on both sides at similar distance from the end of rails. Position crossbar at desired distance from end of table to avoid any interference with table mechanisms and tighten both crossbar positioning knobs to secure.


⚠ WARNING

- This system consists of crossbar with attached rail clamps and dovetailed mounting plate for stabilizer. The rail clamps are designed to fit the rectangular side-rails of standard sized operating and cystoscopy tables.


- I. Rail Clamp Levers
- II. Stabilizer Quick Connect Knob
- III. Rail Clamp
- IV. Stabilizer Mounting Plate
- V. Crossbar
- VI. Crossbar positioning knobs

2. Remove stabilizer from transportation stand. Fit dovetailed mounting plate to mating stabilizer mounting plate on table mount. Secure with stabilizer quick connect knob. Release locking control knob and manually position the fine adjustment mechanism (FAM) in a convenient position to facilitate positioning of fine adjustment controls and attachment of stepper. Tighten locking control knob to secure.

- I. Fine Adjustment Mechanism (FAM)
- II. Locking control knob
- III. Weight adjustment knob
- IV. Dovetail mounting plate

3. Check to see all FAM controls are set near midpoints of their range of travel.

- I. Knob YR - Rotates probe to parallel with prostate centerline
- II. Knob ZR - Levels prostate horizontally on transverse image
- III. Stepper Quick Connect Knob - Secures stepper to FAM
- IV. Stabilizer interface plate
- V. Knob X - Centers probe along midline axis of the prostate and centers the transverse image on grid
- VI. Knob XR - Tilts probe parallel to the posterior surface of the prostate in the sagittal plane
- VII. Knob Z - Used to "zero" the stepping unit at the bladder neck to set the base plane for implantation
- VIII. Knob Y - Adjusts image height on grid to define posterior row for seed placement

NOTE: Knob Y controls vertical movement and has a useful range of about 6cm (3 up and 3 down). Going down it will bottom out as lower plates come into contact and it is shipped in this position. To setup for use, turn Knob Y until lower plates are separated by 2-3cm.

ADJUSTING THE STABILIZER

1. Grasp handles of FAM and loosen locking control knob on stabilizer.

For Micro-Touch® Stabilizer:

Adjust weight by turning weight adjustment knob on stabilizer to achieve desired feel to equipment.

2. Obtain satisfactory transducer position and starting image, tighten locking control knob on stabilizer. After stabilizer is locked into position, locking control knob should not be used until removal from the patient.
3. Use FAM to optimize position of probe.
4. Install sterile drape and grid, position grid platform appropriately and secure with grid rail locking knobs. Refer to Sterile Drape and Template Grid instructions for use.

FUNCTIONAL PREUSE CHECKS

WARNING


- Prior to each use perform the following checks to ensure optimal performance.
- If locking control knob on stabilizer will not hold unit securely in place, please contact CIVCO.

1. Attach the table mount at a suitable clasping point on the table rails to ensure it remains securely in position.
2. Ensure stabilizer moves freely when locking control knob is loosened and secures firmly when locking control knob is tightened.
3. Move each of the knobs on the FAM through its full range and then center each knob in position.

ASSEMBLING TRANSPORTATION STAND


For use with Micro-Touch® Stabilizer LP (Single-Sided Table Mount)

1. Insert threaded stud end of tube through center of rolling base.
2. Slide locating washer over threaded stud and nest into center of hole in rolling base.
3. Turn locking knob onto threaded stud and secure firmly.

- I. Threaded stud
- II. Rolling base
- III. Locating washer
- IV. Locking knob


4. The transportation stand provides a convenient and mobile platform to allow ease of cleaning and handling within the surgical environment. Stand shown with stabilizer and table mount clamped in place.

- I. Stabilizer
- II. Table mount
- III. Transportation stand


SETTING UP MICRO-TOUCH® STABILIZER LP (SINGLE-SIDED TABLE MOUNT)

1. Loosen rail clamp levers on table mount before removing from transportation stand. Slide table mount onto table rail and tighten rail clamp levers to secure. Center the Y adjustment.


- I. Stabilizer Quick Connect Knob
- II. Knob Y
- III. Rail Clamp
- IV. Rail Clamp Levers
- V. Stabilizer Mounting Plate

2. Remove stabilizer from transportation stand. Fit dovetailed mounting plate to mating stabilizer mounting plate on table mount. Secure with stabilizer quick connect knob.

3. Check to see all FAM controls are set near midpoints of their range of travel.

- Knob X - Centers probe along midline axis of the prostate and centers the transverse image on grid
- Knob YR - Rotates probe to parallel with prostate centerline
- Knob ZR - Levels prostate horizontally on transverse image
- Stepper Quick Connect Knob - Secures stepper to FAM
- Stabilizer interface plate
- Knob XR - Tilts probe parallel to the posterior surface of the prostate in the sagittal plane
- Knob Z - Used to "zero" the stepping unit at the bladder neck to set the base plane for implantation
- Knob Y - Adjusts image height on grid to define posterior row for seed placement (Located on the outboard end of the rail clamp.)

NOTE: Due to low clearance of sliding tabletop in relation to fixed components in some tables, a foam mattress pad may be required to allow complete freedom of table movement when the single-sided table mount is in operating position with ultrasound transducer in rectum.

ADJUSTING THE STABILIZER

1. Grasp handles of FAM and loosen locking control knob on stabilizer.
For Micro-Touch® Stabilizer:
Adjust weight by turning weight adjustment knob on stabilizer to achieve desired feel to equipment.
2. Obtain satisfactory transducer position and starting image, tighten locking control knob on stabilizer. After stabilizer is locked into position, locking control knob should not be used until removal from the patient.
3. Use FAM to optimize position of probe.
4. Install sterile drape and grid, position grid platform appropriately and secure with grid rail locking knobs. Refer to Sterile Drape and Template Grid instructions for use.

FUNCTIONAL PREUSE CHECKS

WARNING

- Prior to each use perform the following checks to ensure optimal performance.
- If locking control knob on stabilizer will not hold unit securely in place, please contact CIVCO.

1. Attach the table mount at a suitable clamping point on the table rails to ensure it remains securely in position.
2. Ensure stabilizer moves freely when locking knob is loosened and secures firmly when locking knob is tightened.
3. Move each of the knobs on the FAM through its full range and then center each knob in position.
4. Perform table clearance check. Before indexing the table, make sure stabilizer and table mount do not interfere.

REPROCESSING **WARNING**

- *Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.*
- *See your system's user guide for reprocessing transducer between uses.*
- *These procedures for reprocessing have been validated for effectiveness and compatibility. Equipment could be damaged or cross-contaminated due to improper reprocessing.*
- *Do not place in a mechanical washer.*
- *Do not use ethylene oxide or autoclave to sterilize.*

1. Disinfect all surfaces with common germicidal or antiseptic solution, such as alcohol or hydrogen peroxide.
2. Wipe dry with lint-free cloth or allow to air dry prior to reassembly and use.

STORAGE AND DISPOSAL **WARNING**

- *Dispose of single-use components as infectious waste. Reprocess reusable components after each use according to instructions for use.*

- NOTE:**
- Device should be stored to keep all components together and secure.
 - For questions or to order additional CIVCO products, please call +1 319-248-6757 or 1-800-445-6741 or visit www.CIVCO.com.
 - All product to be returned must be in its original packaging. Contact CIVCO for further instructions as needed.

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**Australian
Sponsor Address**

SCANMEDICS
Unit 6, 15-17 Gibbes St.
St. Chatswood, NSW 2067
Australia

CIVCO Medical Instruments Co., Inc.
102 First Street South, Kalona, IA 52247-9589
319.248.6757, 800.445.6741
Fax: 319.248.6660, 877.329.2482

COPYRIGHT © 2020 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ACCUCARE AND MICRO-TOUCH ARE TRADEMARKS OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW.

WWW.CIVCO.COM

2020-07-27 043-659-3A

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MPS GmbH